

KX-79-1 MAGNETISM HEATING MIXER

How to use

Using this instrument before, first, check whether the attached accessories are complete, and put the beaker and put the mixer in it. Then, plug the power device: connected to the power and open the power switch, indicator light that is to work. Just turn the heating switch to the heated stirring, with no mixing, when not heated. There is slow speed too fast (Files do not allow fast start, avoid the mixer to jump). Cut off the power when not working. To ensure safety, plug in the ground, and the equipment should be kept clean and dry, non-solution into the machine to avoid damage to parts.

Introductions

This instrument is widely used in the laboratories of major colleges environmental, protections, scientific research, public health, epidemic prevention, petroleum, metallurgy, chemical industry, medical and so on. The laboratories must have the ideal tool to greatly enhance the efficiency of the experimental staff.

Characteristics

The machine uses permanent magnet DC motor, no-load can be operated more than 10,000 hours continuously. The instrument performance is good, no noise, no vibration, mixing results are obvious, for laboratory staff, providing a more ideal tool. This machine is cold Rolled Plate touch welding forming. The outer surface of high temperature spray enhances the corrosion resistance.

Capabilities

1. Power supply: 110V 60HZ
2. Motor speed: starting:2400r/min (stepless adjustable)
3. Motor power: 25W
4. Hot-wire power: 300W
5. Heating: 200CC solutions at room temperature of 20 °C, heat 8 minute or so, usually ranging from 90 °C

Note:

1. Found the mixer sub-mixing or jumping, please cut off the power to check whether the beaker is stable or location. In the meantime, please test whether the now-used voltage between 110V \pm 10V, otherwise there will be above.
2. Heating time in general should not be too long, intermittent usage can extend the service life.
3. Medium-speed operation can continuously work for 8 hours, high-speed operation can continuously work 4 hours.
4. The suite of machines is equipped with a set of clamps and pillars a fluorine mixed, a power cord (plug). To ensure credibility, those who buy the products from our factory find quality problems caused by failures themselves within one year from the production date. We promise to repair or exchange it for you.